



WET ABRASION SCRUB TESTER

CONFIRMS TO ASTM D-3450,2486,4213,4828; DIN-53778; BS7719,3900 F-17; ISO11998;

One of the most basic equipment to test the abrasion resistance / durability / gloss retention / endurance / cleaning properties etc., Of Detergents, Paints, varnishes, inks, plastics, Floor coverings, carpets, upholstery materials. Ceramics - glazed or unglazed, Tiles, Furnitures, Leathers, Garments, Wallpapers, etc.,



Features :

- ✎ 4 digit programmable pre-set counter with LED display and battery backup for unattended testing
- ✎ 2 Scrubbing heads and tracks allow carrying out of Simultaneous Comparative tests.
- ✎ Self Cooling TEFC motor for trouble free performance
- ✎ 3mm Glass plate on a Stainless Steel top tray make a non-corrosive test surface
- ✎ Adjustable Flow rate of test solution by a manually controlled Dripper / Programmable and precise peristaltic Pump

Models :

WST-B - BASIC INSTRUMENT WITH MANUAL SOLUTION DISCHARGE ARRANGEMENT.

WST-D - WITH BATTERY BACKUP AND PERISTALTIC PUMP FOR PRECISE DISCHARGE OF SOLUTION.

Standard Accessories : 2 nos. Brushes, Carriage Assembly and weights to test as per any one of the IS test methods.

Optional Accessories : Different types of Brushes, Suitable Carriage Assemblies and weights to test as per any of the above mentioned national and international standards.

Due to regular development, the design and specification of the instrument is subject to change without prior intimation.



Works Address : D-428, 10th Main, Peenya II Stage, Peenya Industrial Area, Bangalore-560 058.
Telefax:091-080-23364115/23349239
E@mail: