



STIRRERS

MAGNETIC STIRRERS

CODE:22001

Stainless steel top housing series of magnetic stirrers are compact, with PMDC motor, lightweight, accurate speed control up to 1200Rpm, Ideal for laboratories, schools etc., Easy to handle.

MODEL **CAPACITY**

1 ML	1 ltr.
2 ML	2 ltr.
5 ML	5 ltr.

CODE:22001-01
22001-02
22001-03



Complete with stir paddle and manual.

‘REMI’ MAGNETIC STIRRERS WITH HOTPLATE

These are similar to ML series but with additional Stainless steel top hotplate. They are compact, with PMDC motors, lightweight, accurate speed control up to 1200Rpm, Ideal for laboratories, schools etc., Easy to handle.

MODEL **CAPACITY** **HEATING**

1 MLH	1 ltr.	150W	22001-04
2 MLH	2 ltr.	300W	22001-05
5 MLH	5 ltr.	500W	22001-06



Complete with stir paddle and manual.

‘REMI’ MAGNETIC STIRRERS - PLUS SERIES

- * Speed 100 to 1500Rpm,
- * They are compact, lightweight, and accurate.
- * Porcelain Enamel Hot plate top for upto 340°C
- * Ceramic Hot plate model for up to 550°C
- * Backlit LCD display of set and actual parameters
- * External temperature probe (PT-1000) for control of fluid temperature
- * Scratch Proof and chemical resistant
- * Brushless motor, with Protection class to DIN En60529.

MODEL **CAPACITY**

5 ML Plus	5 ltr.	135mm dia.	22001-07
10 ML Plus	10 ltr.	184mm x 184mm	22001-08
5 MLH Plus	5 ltr.	135mm dia.	22001-09
10 MLH Plus	10 ltr.	184mmx 184mm	22001-10



Complete with stir paddle and manual.

‘ANTECH’ MAGNETIC STIRRERS WITH HOTPLATE

- * Speed Up to 1500 Rpm,
- * They are compact, lightweight, and accurate.
- * Ceramic coated Hot plate model for up to 280°C
- * LED display of set and actual parameters
- * External temperature probe (PT-1000) for control of fluid temperature
- * Scratch Proof and chemical resistant
- * Easy to handle. Specially designed for corrosive environment.

MODEL **CAPACITY**

AN-MSH-680	2 ltr.	22001-11
------------	--------	----------

Complete with stir paddle and manual.